

Round CollectionTM

by EUROSILICONE®



Eurosilicone® is a leading global manufacturer of innovative quality medical devices.




Overview

- Currently features 12 different configurations with 233 sizes.
- ✓ Different projections: Low, medium, medium extra-high, high, extra-high.
- ✓ Shell surface:
 Smooth or Cristalline™ Microtextured.
- ✓ Silicone gel Paragel™: Soft Cohesive or Natural Cohesive.
- ✓ ParagelGard™ 360° Barrier Layer: Eurosilicone's barrier layer technology minimising qel diffusion.

Natural Cohesive Design Microtextured Shell



Soft Cohesive Design Microtextured Shell



Soft Cohesive Design Smooth Shell





- Duteille et al. Eight-Year Safety Data for Round and Anatomical Silicone Gel Breast Implants. Aesthetic Surgery Journal sjx117. 2017, 1–11. DOI: https://doi.org/10.1093/asj/sjx117.
- H. Hoflehner, S. Spendei, D. Lumenta, L. P. Kamolz. Das Österreichische Register für Silikonimplantate. Dermatologie & Plastische Chirurgie. 4/2014; 26-32 [The Austrian Silicone Implants Register. Dermatology & Plastic Surgery. 4/2014; 26-32].
- ³ Internal data, % calculated on the basis of the number of complaints received compared to the number of units sold during the same period
- * Actual rates per implant, all cohorts were pooled in this analysis.
- ** Reoperation rates for complications deriving specifically from the product.

Why choose the Round Collection™

The Round Collection provides a wide range of products, including smooth and microtextured surfaces, two types of gel cohesivity, and four different projections – from low to extra-high.

Paragel™

Paragel is the high cohesive, medical grade gel, used in all Eurosilicone implants. It provides controlled gel distribution and is available in both Soft and Natural cohesivities; giving more options to surgeons and patients.

ParagelGard[™] 360° Barrier Layer

ParagelGard™ 360° is Eurosilicone's protective barrier, which has specific technology designed to minimise gel diffusion. The barrier is also present at the product patch, completely covering the implant.

Cristalline[™] Texture

The Round Collection uses the microtexture silicone shell surface, which is designed to offer tissue adhesion, minimising rotation and capsular contracture.





Proven Safety through Clinical data.

Eurosilicone breast implants have been used by physicians around the world for almost 30 years. During this time, long term performance data has been collected via clinical studies, as well as internal safety data, which demonstrates both high quality and safety of our product range.



Ongoing Post Market Clinical Study.

Capsular Contracture Rate¹

Rupture
Rate¹

**The state of the state of

Reoperation¹ **

Austrian Breast Implant Registry

Number of implants in Austrian Registry² 2004-2014: 17540.

Capsular Contracture²
Total number of implants: 2058.

Confirmed Rupture²
Total number of cases: 444.

2629
Other: 15%

2350
Allergan: 13%

6618
Mentor: 38%

5943
Eurosilicone: 34%

731
Other: 35.5%

749
Allergan: 36.4%

317
Mentor: 15.4%

261
Eurosilicone: 12.7%

96
Allergan: 21.6%

90
Mentor: 20.3%

Internal Safety Data Report

Total Capsular Contracture Rate³

0.45%

Total Rupture Rate³

0.2%

Relax, you're covered...

Introducing the GCA Comfort Guarantee

As a GCA breast implant, the Round Collection offers the GCA Comfort Guarantee, a lifetime guarantee against rupture and capsular contracture.





For more information about our products please visit **www.gcaesthetics.com** or email ComfortGuarantee@gcaesthetics.com.







